



Be Right™



TEST KIT, 5-EP mg/L HARDNESS

Produktnr:

145401

Farligt

SEK Pris:

Kontakta oss

Tillgänglig



Total Hardness drop count titration/EDTA test kit with range of 20 - 400 mg/L; contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specificationer

Antal tester :	100
Kapslingstyp:	B
Mätområde:	20 - 400 mg/L CaCO ₃
Metodnamn:	Drop count titration/ EDTA-powder pillows
Parameter:	Hardness, total - As CaCO ₃
Plattform:	Drop Count Titration
Smallest increment:	20
Test-kit modell:	5-EP MG-L