



## Test kit, hypochlorite, model CN-HRDT

**Produktnr:**

**2687100**

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Hypochlorite Kit, Model CN-HRDT with Digital Titrator

Utilizing the Hach Digital Titrator, this kit is used to determine hypochlorite, high range (Bleach) solution concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 5 - 15% as Cl<sub>2</sub>: 50 - 150 g/L as Cl<sub>2</sub>. Contains 100 tests.

**Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration**

**Rugged titration device is ideal for field use, with no glass burets required**

**Designed for on-site analysis of bulk hypochlorite deliveries or bulk storage solutions at drinking water or chloramination facilities**

**Verify concentration of hypochlorite solutions at time of delivery**

**Kit contains all required reagents, solutions, and glassware**

### Specifikationer

Antal tester :	100
Kapslingstyp:	G
Mätområde:	5 - 15 % = 50-150 g/L Cl <sub>2</sub>
Metodnamn:	Digital Titrator/ Thiosulfate
Metodnummer:	10100
Parameter:	Hypochlorite, high range (Bleach) - As Cl <sub>2</sub>
Plattform:	Digital Titrator
Smallest increment:	0.05%
Test-kit modell:	CN-HRDT