



Fast volym, mikropipett, 20 µL

Produktnr:	2946003
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Micropipette features a two-stage plunger dispenser, which ensures high accuracy and precision in pipetting techniques and facilitates last drop dispensing. Unique tip cone accepts both micro and ultra micro pipet tips. Each micropipette is fully autoclavable and color-coded for easy identification. Includes one pipet tip. Additional pipet tips sold separately.

Specifikationer

Kapacitet :	20 µL
Noggrannhet:	0,5 %
Volym/förpackningsstorlek:	1 each